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CHOICE Dr. **DC-001**

CHOICE Dr.

Thank you for choosing and supporting CHOICE Dr.

We hope that you are satisfied with our product and service, and that we have the opportunity to serve and accompany you and your family.

If you have any feedback or reviews regarding the quality of our product and service, please contact our hotline at **1800 8379** so that we can receive and promptly address them.

In addition, please scan the **Zalo OA** code to activate the electronic warranty, receive additional gifts, and enjoy promotion and **24/7** support from CHOICE Dr.

Thank you once again and we are happy to have take care of you and your family's health!

Choice Dr

CONTENT

NOTE

OPERATING INSTRUCTIONS

SAFETY INSTRUCTION

ADULT ONLY

SPECIFICATIONS

MASSAGE HEADS

1. Round

It is the acceptance and push between the type of propagation and diffusion of the output of the power generated by the spherical design. Provides even and soft pressure. For massage and relaxation of the waist, abdomen, hips and legs. It is recommended for lightweight training lovers for dancing, yoga, aerobics, etc., besides muscle relaxation before and after exercise.

2. Mushroom-shaped

It features a small and simple design. In addition to the flat surface, small curves have been added to fit the mechanism of the body shape for a better feel. For use on the internal and external muscles of the abdomen, chest and legs of the human body.

3. U-shaped

With even force between the two points, the stimulation is strong and effective for simple acupuncture. A reasonable u-width dimension perfectly bypasses the human spine and cervical vertebrae for a better experience. For use on medical, cervical and bilateral muscles. (Should not be used for direct impact on the cervical spine and spinal column as this may cause physical injury)

4. Tapered

It has a simple shape but is packed with impactful experiences and applies pressure directly to the depths. For use on the meridians, soles of the feet, palms of the hands, Achilles tendons and rhomboids of the back.

OPERATING INSTRUCTIONS

1. Power on/off button (a long press to turn on/ a long press to turn off)
2. Press the on/off button to change gears
3. Removable and replaceable vibrating head. To remove the vibrating head, simply pull it out by force
4. Please choose a 5V charging device for charging.
5. Please charge for 3 hours for the first time, and connect the SV adapter with the data cable to charge.
6. button model: the 6 lights flash alternately when charging, the 6 lights are always on when fully charged, and the lights indicate the gear when working.
7. Touch model: the digital tube shows the power level when charging, and shows the gear level when working.

SAFETY INSTRUCTION

To minimise the risk of electric shock, fire and personal injury, it is essential that you use this product by the following instructions.

ADULT ONLY

- This product should only be used on a dry, clean skin surface and should not be used through clothing. press and move gently over the skin. The duration of each part is approximately 60 seconds.
- Only use the device on the soft tissues of the body without causing painful discomfort. Do not apply it to any hard or bony parts of the head and body.
- Only use the vibrating head that produces the best effect.

- Bruises can be caused by either moderate pressure or normal pressure. Please always check the area being massaged.
- As soon as there is pain or signs of discomfort, discontinue use.
- Finger hair or other body parts should be kept away from the shaft and rear of the vibrating head as they may be crushed.
- Do not put any objects in the ventilation holes of the appliance. Do not immerse it in water or let water enter the ventilation holes of the appliance. Do not drop the appliance or use it improperly. Use only 5 VDC to charge the appliance charger.
- Carefully check the appliance and battery before each use. Do not alter or modify the appliance in any way.
- Never operate or charge the appliance in an unattended manner.

SPECIFICATIONS

Model: DC-001	Battery: 1200 mAh
Voltage: AC240V/50HZ	Weight: 0.5 kg
Power: 30W - 0.8A	Product size: 13*5*14cm