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EN

IMPORTANT SAFETY INSTRUCTION

Danger: In order to reduce the risk of electric shock, fire and personal injury, this product must be used in accordance with the following instructions:

- For adults only it can only be used on dry, clean skin surfaces of the body, and cannot be used through clothes. Press and move it gently on the skin. The duration of each part is about 60 seconds.
- Use this device only on soft tissues of the body without causing pain or discomfort. Do not use it on the head or any hard or boned part of the body.
- Only use the vibrating head that can produce the best results.
- Bruise may be caused whether under moderate pressure or normal pressure. Please always check the area that receives physiotherapy.
- Stop using it immediately once there is any sign of pain or discomfort.
- Keep fingers, hair or other parts of the body away from the shaft and back of the vibrating head, because they might be pinched.
- Do not place any object in the vent of the device.
- Do not immerse it in water or allow water to enter the vent of the device.
- Do not drop the device or use it improperly.
- Only charge the device with the supplied charger. Do not use another charger to charge it.
- Carefully check the device and battery before each use.
- Do not change or modify the device in any way.
- Never operate or charge the device in an unattended way.

INDICATIONS AND CONTRAINDICATIONS

The deep penetration of the device can reduce pain from joint sprains, bursitis, muscle cramps, neuralgia, and other musculoskeletal disorders in a short period of time. Almost all soft tissue diseases are associated with trauma, excessive muscle use or abnormal posture, and most of the scar tissues involved in these diseases can be treated with this device. Avoid excessive stress and prolonged physiotherapy when using the device. Also avoid sensitive areas such as the head, face, and cervical vertebrae. Stop physiotherapy if there is inflammation, swelling, or increased pain. The ongoing research is about the impact of vibration training on specific medical conditions. This is very likely to reduce the types of contraindications mentioned above, and practical experience has shown that in many cases, it is also a physiotherapy recommendation to combine vibration training into physiotherapy programs. This must be done by, or in the company of, a doctor, expert, or physical therapist.

1

Indications	Contraindications
Pain and cramps due to muscle injury, sprains, strains. Help the flow of edema fluid in the swollen area. Relax thickened connective tissues and fascias. Reduce the accumulation of lactic acid in muscles. Increase joint mobility. Eliminate muscle fatigue.	Aneurysm, bleeding, use of blood thinner. Heart disease, with a pacemaker or defibrillator. Pregnancy, cancer. Within 90 days of intra-articular fixation, the device should not be used within 3 inches of internal fixation. Sensitive area: head, face, cervical vertebrae spine Positions close to the bones with less muscle coverage: tibia, the back of the foot, the back of the hand.

CONFIGURATION AND SPECIFICATIONS

2

OPERATE INSTRUCTION

Using method

- Insert the suitable vibration head into the massage gun.
 - U-shaped massage head: suitable for cervical spine, neck, and Achilles tendon.
 - Ball massage head: suitable for large muscle groups such as the spine and thighs.
 - Bullet massage head: applicable to the meridians, joints, palms, and soles of muscle in various parts.
 - Flat massage head: applicable to all parts of the muscle to relax, shape, and beautify body lines.
- Press the Power button to start massage. The massage head will start to vibrate.
- Press the power button again to choose the suitable massage speed for your body. The red light will increase after each pressing time indicating the current massage speed.
- Bring the massage head to the position need to be massaged, then move the massage gun slowly up and down or left and right.
- Recommend using time: 60s for each part, maximum 15 mins/ day.
- After finishing the massage, press the power switch until the massage gun is off.

Operate program for lower limbs

Flat massage head vertically downward

- Hamstrings include biceps femoris, semitendinosus, semimembranosus: 20-30s x 3 times.
- Quadriceps femoris include vastus lateralis, vastus medialis, vastus intermedius and rectus femoris: 30-45s x 2 times.
- Sartorius: 30-45s.
- Tensor fasciae latae, tibialis anterior, peroneus longus, peroneus brevis: 15-20s x 2 times.

Flat massage head outward 45°

- Gracilis, adductor longus, adductor brevis: 10-15s x 2 times.
- Triceps surae, gastrocnemius and soleus: 20-30s x 3 times.

Bullet massage head vertically downward

- Popliteus: 15-30s.
- Iliac bone edge – Iliac bone anterior: 30s – 1 min.

U-shaped massage head vertically downward

- Popliteus: 15-30s
- Triceps surae, gastrocnemius and soleus: 20-30s x 3 times.

Note: The movement direction of massage head is from up to down.

3

MAINTENANCE, CLEANING, AND STORAGE

Clean the device with a slightly wet towel, and then wipe it dry with a soft cloth for storage, or put it into the supplied carton when you go traveling.

WARNING AND REMINDERS

Do not immerse it into water. Keep it away from liquids or heat sources. Keep the vent free from dust and scraps.

Do not remove the screws or try to disassemble the device.

Do not continuously use it for more than 1 hour. Please stop the device for 30 minutes and then use it again.

Unplug the charger from the battery after charging or before using.

It is only for adults. If you have injuries, do not use this device. Please consult with your doctor before using this massage gun.

CHARGING

Before using the massage gun for the first time, please charge it fully for 5-8 hours. And then, for subsequent charges, you only need to charge it for about 4 hours.

During charging, connect the type C end of the supplied adapter to the charging port of the battery, and then plug the cable at the USB side onto the wall outlet

When the battery is connected, the LED on the device flashes, in dicating that the battery is being charged

When all LEDs stop flashing and are on, it indicates that the battery is fully charged.

The battery can be charged at any time and at any level.

The average operating time of the device is two hours, depending on the speed and the applied pressure during use.

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