

### Introduction

Thank you for purchasing electronic scale. This instruction manual contains important information about safety, usage, and disposal. Use the product as described and keep this manual for future reference. If you sell this scale or pass it on, also give this manual to the new owner.

### Safety Instructions

- Do not immerse the scale in water or other liquids.
- Do not place the scale on carpets or soft surfaces. Only place on solid and flat surfaces.
- Do not step on the scale with wet feet and do not use the scale if the platform is moist. Moisture will cause the scale's surface to become slippery.
- Do not jump on the scale as this may damage the sensors.
- Do not step on the edge of the scale while getting on and off, as it may cause the scale to tip over.
- When using the scale for the first time, or when changing its location, it should be turned on, off, and on again. Its surface should also be free of any objects. Failure to do this may result in inaccurate measurements.
- Always keep batteries out of the reach of infants and young children to prevent them from being swallowed. If swallowed, take immediate action and contact medical authorities for help.
- Clean the scale with a slightly damp cloth. Do not use any chemical/abrasive cleaning agents.
- If you do not intend to use the scale for a long period of time, remove the batteries and store the scale in a clean, dry place away from direct sunlight.

### Inserting Batteries

To insert the batteries, open the battery compartment on the rear of the scale and remove the compartment cover. Insert 3 AAA batteries, paying attention to the correct polarity.

### Select Unit Before Use

1. Wait until "0.00kg" appears.



2. Press the UNIT button on the rear of the scale to select the desired unit (kg, lb or st).



Note: After selecting the unit, place the scale on a solid and flat surface and start the measurement as described below. The scale will shut off automatically after 8 seconds of inactivity.

### Re-calibrating the Scale

Re-calibrate the scale before weighing if the scale has been physically moved or the batteries have been changed.


Directions: Step on the scale, wait for the display to turn on, then step off. Wait until the scale automatically drops to 0.0kg/0.0lb and step back on the scale and begin weighing.

1

2

### Measuring

1. Place the scale on a solid and flat surface.
2. Stand on the scale to turn it on.



3. Stand still and the display will blink three times presenting the confirmed weight measurement. The scale will then turn off after 10 seconds.



Note: If you want to measure again, please step off completely and then step back on the scale.

### Tips for Tracking Accurate Results

1. Weigh yourself in the morning (either daily or weekly) to get your weight before you start the day.
2. Make sure to re-calibrate the scale before weighing if the scale has been physically moved or the batteries have been changed.
3. To receive the most accurate results, please remove all articles of clothing and jewelry.
4. Only use the scale to measure your progress. Ideally, users should consult with their doctors in settling long term goals.
5. Make sure to test yourself at the same time of the day under the same conditions.

### Maintenance and Care


- Simply use a soft damp cloth or sponge to wipe any dust on the scale.
- DO NOT get any water or liquids into the electrical systems inside.
- For stains or residue, use a cloth with rubbing alcohol to remove it.
- Note: This scale is not waterproof, DO NOT submerge in water.

### Indication Messages

- "Lo" is shown when battery power is low. The scale will then turn off after a few seconds. Please replace with new batteries.



- "Err" is shown when the scale is overloaded. Please remove the weight on the scale immediately to protect the sensors.



### Product Specification

Scale Size: 280×280×25mm  
 Units: kg/lb/st  
 Capacity: 0.1kg-180kg  
 Division: 0.05kg/0.1lb  
 Power: 3×1.5V AAA

Starting Weight: ≥5kg  
 Working temperature: 5~40°C  
 Working Humidity: ≤90%RH  
 Storage Temperature: -10°C~60°C

### Device Requirements

iOS devices require iPhone 4S, iPad and above.  
 Android device system requires Android 4.3 and above.

### Precision Accuracy

50kg: ±100g 100kg: ±200g 150kg: ±300g

3

4